

December 18-21, 2020

(v13) HEU: DTX/CRM197 Optimization "#2"

Overnight rest, 2 hours stimulation, 18 hours incubation;  
+/- Costim, 5, 10, 20  $\mu$ g doses.

Costs w/ reg 60% line in Million calls (as volume)

Specimen	Status	Location		Date	Tasks	Volume	Ly	Ly+Mon	Total	
ND050	Adult T1/T2	-130°C	Rested	10/21/20		2.5	2.33	3.60 E4	5.75	2.875
ND050	Adult T1	-130°C	Fresh	3m	LSIE4/2.44	-1 ml (warm)	4.89	6.35 E4		

Adult 1.0E+6 cells	Ex vivo	Extracellular * 3.5		Intracellular * 3.5	
V450 BV421	CD69	1	3.5	CD69	1
V525 BV510	Aqua	1:500		Aqua	-
V670 BV650	CD4	1.8	6.3	CD4	0.5
B530 FITC	CD8	1.5	5.25	CD8	0.5
B710	TNFA	-		TNFA	3.0
PerCPeF710					
Y590 PE	CD62L	2	7.0	CD62L	-
Y780 PE Cy7	IFNg	-		IFNg <sub>x2</sub>	0.8
R670 APC	CD3	1.5	5.25	CD3	0.5
R780 APC- Fire750	CD45RA	1.5	5.25	CD45RA	-
PBS		10.2	39.2		13.7

☐ running low antibody, waiting delivery

5ml - 10  $\mu$ l  $\checkmark$  0 Aqua

2/10 @ 9:05 am

FB 5 pm 9:23 am

Surface @ 9.35am

G/H/P @ 10:09 am

Perm wash @ 10:34 am

Intra @ 10:53am

Done @ 11: am

Thawed @ 5:35pm  $\rightarrow$  10 ml media/2 wells.  
Processing started @ 11:00 AM (17 1/2 hours rest)

add wash @ 11:18 → very few cells left mostly  
adhered - Give 10 minutes on ice today.

In centrifuge @ 11:30 am (still the cell) random  
cell no more than 3 per field (view)


Stack coat @ 11:54 a.m.



13c/ME ad  $\rightarrow$  Virtually recta. 4 lymphocytes


Hot - 100


180C/µm ND050 → [Thaw another for SD's:]

Nerve, TI Csm R7

$+ 2k_{\text{eff}} \text{ catim}$ 


$+ 1k_{\text{eff}} \text{ catim}$ 




Inoculation start @ 12:49 pm + 20h = 9 AM

Let's use a scaled golgi block/plug approach this time?

1 + 1 + 19 ~ 10  $\mu$ l ~ 1M / 1ml  
in 500  $\mu$ l mother  $\leftarrow$   
vs 1M cells  $\leftarrow$

Let's stick w/ what works I guess...  
everyone gets 10pt for now, repeat  
w/ 5pt in January,

[2.49 pm] Add Goli Bldl.

13:02 Thaw new ND050 for sc's etc.

13:20 Ine Cell count

Bit more debray field on FSCX 55r.

4.53 M Total 1.51 M/ml

$$300 \cancel{\text{K}} \cdot 1.51 \text{ M} = 0.198$$

14.04 sc's inc start

2:20 put leftover RNA in freezer

2:39 sc's done,

Cell in for bio med @ 14:49 pm  
So back by 8:49 AM

500  $\mu$ l media  $\sim$  1 million cells  $\sim$   $\frac{\text{Concentration}}{2M/ml}$

$\Delta$  Golgi Conc

10  $\mu$ g/TT/ml  $\cdot$  million cells  $\rightarrow$  5  $\mu$ g in 500  $\mu$ l  $\sim$  1.785  $\mu$ l  
 10  $\mu$ g/Crm199/ml  $\cdot$  million cells  $\rightarrow$  5  $\mu$ g in 500  $\mu$ l  $\sim$  0.5  $\mu$ l of stock.

~~10  $\mu$ g~~

Costim 5 + 5 + 42  $\mu$ l = 52  $\mu$ l  $\rightarrow$  0.5  $\mu$ g/ml  
 so 12.5  $\mu$ l for a 500  $\mu$ l volume

Costim	Name	TT	Crm199	TT
k				
2x				

Golgi Plug/Golgi stop

1 + 1 + 19  $\mu$ l = 21  $\mu$ l for 1M in 1ml?  $\leftarrow$  Previous focus # cells  
 5  $\mu$ l for 1M in 500  $\mu$ l  
 which factor matters most

3M single colors

5  $\mu$ l - 50  $\mu$ l  
 1:10  
 10  $\mu$ l = 0.5  
 5  $\mu$ l = 0.25  $\times$  2(500  $\mu$ l) = 0.5  $\mu$ g/ml effective 2x

18 ml 8% FBS - PBS  
 10%

1.8 ml  
 900  $\mu$ l - into 17.1  $\mu$ l through streamer.  
 10 ml 1 ml mb 19 ml