

November 10th, 2022

ILT Panel Final Checks #1

1:100

Specimen	Status	Location	Conc	Date	Tasks	Volume	Ly	Ly+Mon	Total	2E+6	.2E+6
ND006	adult w/2 adler 200 adler		15E+6	Mar 11-2022		3	2.32	3.23	6.92	862 μ l	86.2 μ l
ND006	adult w/2 adler		15E+6			2	3.72	4.88	7.44	537 μ l	53.7 μ l
									14.36		

stain count @ 12:38 pm
Count on beaver @ 12:51 pm

2 μ l PMA + 500 μ l + 500 μ l

○ ○ ○ ○ (2M)

→ 8M (6M/30)

Spiked w/o → from actual sample

CD107A + IFNg + TNFa + CD69 + unstained / → 1 μ l PMA + 500 μ l (131)

CD62L unstained

→ 200K cells

4 million act + 1.5 tubes (act controls)

4 million act + 1.5 tubes ch controls

PMA/ctrl in @ 1:37, CD107a @ 6 μ l w/ samples

06:02 pm decided adjusted ab stain amounts to test
06:03 pm ab spin down
Abs prepared by 6:54 pm

5
x.8
4.0 → 2.0 μ l L/O NIR in 5 ml

6 * 4 + 2 cont = 36 ml 5% FBS

30 * 2 = 60 ml sc's FBS

46 ml → 100 ml 5 ml

Rb lysis

300 μ l * 6 = 1.8

* 30 = 9.0

11000
1100
9000

6 * 6 = 36 ml Perm wash

30 * 2 = 60 ml

46 ml Perm wash

86 * 100 μ l = 60 μ l

30 * 50 μ l 1500 μ l

2.100 ml 0.4% PFA-RBS

7:25 pm reagents prepped

intra done
@ 23:20 pm
Cold Rb lysis
10:42 to 11:50
@ 22:52 pm

Simple Fix Perm @

11:01 → 11:21

Sc's only got

1 2 ml Perm wash

Δ size correctly?

and sufficient for

intracellular sc's?

1st sample perm wash

1 ml @ 23:21 pm

2nd wash @ 23:32 pm

Intra @ 23:42 →

00:22

Done @ 12:30 pm

no direct lyse for
TNFa/IFNg?

Spin down @ 7:54 pm (20 min)

L/O @ 20:09 → 20:24 pm

2 μ l 1:10 hcdid

5 μ l 1:10 hMRI

Tets @ 20:43 →

21:23

Cold sc's @ 9:00 pm → 9:30

Hot sc's 9:06 pm → 9:36

CD69 .5 @ 4°C @ .5 @ intra

2 combi w/ Rb lysed & fixed

21:29-30 Rb lysis → 33 → 36

Hot Rb lysis @ 38

Hot Samples @ 9:41 → 10:11 pm

Sc's Fix Perm @ 9:57 → 10:07 → 10:17

Cold Samples @ 10:20 → 10:50

Perm wash span 10:26 pm for sc's.
Intra single colors @ 22:46 pm

3
463
x 1.5
2315
2315
4630
69.4.5
-463.0
261.5
238.5

November 10th, 2022

ILT Panel Final Checks #1

#	Filter	Single color (ul)	Ref ctrl	Unmixing ctrl	Fluorochrome	Marker	Clone	Vial lot #	During	4	L/D	Tetramer	37°C (RT) for	4	Surface staining	4	Lysing	Intra mix	3
1	UV2				BUV995	CD82L	(DRG-56)		stain		15 min	(RT) for 40 min	30min		40C for 30min		then	40min intra mix	
2	UV7				AF	AF-UV6													
3	UV9				BUV496	CD8	(BPA-18)												
4	UV10				BUV563	CD69	(H450)												
5	UV11				BUV615	CD44	(1G1)												
6	UV14				BUV661	VE2	(B6)												
7	UV16				BUV737	CD63	(1C6)												
8	UV1				BUV805	CD4	(SK3)												
9	UV3				BUV421	CD127	(A01905)												
10	UV5				Pacific Blue	CD14	(M5C2)												
11	UV7				Pacific Blue	CD19	(H1B19)												
12	UV10				BUV480	CD161	(R6A31)												
13	UV11				BUV510	CD45RA	(H100)												
14	UV13				BUV605	CD56													
15	UV14				BUV711	CD72													
16	UV15				BUV750	IFMγ	(B27)												
17	B2				BUV786	CCR6	(11A9)												
18	B4				AlexaFluor 488	hCD1d													
19	B6				Speck blue 550	CD3	(SK7)												
20	B8				PE	NKG2D													
21	B10				PE-CF594	CD26	(M-A281)												
22	B13				PE-CF5	CD25	(M-A251)												
23	R1				PerCP-Cy5.5	TM6a	(MAA811)												
24	R2				APC	CD16													
25	R4				AlexaFluor 647	hMIR1													
26	R6				APC-R700	CD107a	(H4A3)												
27	R7				Zombie NIR	U/D													
28	R8				APC/Fire 750	CD27	(O223)												
29	R8				APC/Fire 810	CD38	(H172)												
And UNSTAINED CONTROLS !!!										Antibody Total									
										R10 Media									
										Pipette draw volume									
										Pipette draw volume /sample									
										19.5									
										6.0									
										24									
										Antibody Total									
										Brilliant Stain									
										Pipette draw volume /sample									
										61									
										1.6									
										6.4									
										200									
										50.0									
										66.4									
										69.6									
										5.75									
										50									
										150									
										54.75									

Simplified Protocol

1. Aliquot cells 3.0x10⁶ Cells/tube
Raise R10 + cell volume to 500 ul
Add 500 ul of PMA-ionomycin/Ctrl
Incubate 6 hours @37C

Wash with 2 ml PBS, 1400 rpm, 6min

800 ul Live/Dead @RT for 15min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Add Tetramers, incubate @RT for 40min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Add hot Surface mix, incubate @37C for 30 min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Add cold Surface mix, incubate @ 4C for 30min

Add 300-500ul 1x RBC lysis for 3 minutes

Wash 2 ml 5% PBS-FBS 1500 rpm, 6min

300 ul Cytoflu/pern, incubate @ 4C for 20min

Wash twice 1 ml Pern/Wash 1500rpm, 6min

Add Intra Staining, incubate @ RT for 40min

Wash once w/ 2ml Pern Wash

Resuspend 100 ul 0.4% PFA-PBS, store 4C

alkaline sbs for CD15
and AF488/647 using Beckin

CD80
CD86
CD137
CD137a
CD137b
CD137c
CD137d
CD137e
CD137f
CD137g
CD137h
CD137i
CD137j
CD137k
CD137l
CD137m
CD137n
CD137o
CD137p
CD137q
CD137r
CD137s
CD137t
CD137u
CD137v
CD137w
CD137x
CD137y
CD137z

November 10th, 2022

ILT Panel Final Checks #1

12:32pm start


ctrl unstained ~ 1400 events.

PMA unstained ~ 900 events (no monocytes for all effective purposes)

Sels @ 400-
470 events/sec

1:07pm
→ 2:10pm

- 200K ~ 50K



4	UV7	V7	<u>B3</u>	1	4
(3)	↓ UV7	↑ V7	= B3	2	
(2)	↓ UV7	= V7	= B3		3
(1)	↓ UV7	↓ V7	↓ B3		

PMA vs ctrl monocytes similar UV7, V7-A, B3-A

(residual)

(little lower in PMA)

"not dying" vs ctrl monocytes "not dying" slightly higher, ghosts
deads share ctrl monocyte AF signature, tail out right to ghost location,
These ghost represent the most fully AF (edge/extreme pericarp)

unique to PMA well

PMA Lymphocytes very similar in UV as controls, but from V2 → V10 have an increase

PMA edge (lesser extent ctrl)

↑ UV7 ↑ V6/V8 ↑ B3

ctrl unique

↓ UV7 ↑ V5/6

November 10th, 2022

ILT Panel Final Checks #1

Checklist for ILT Panel Final Checks #1

1. Check for missing components

2. Check for correct wiring

3. Check for correct labeling



4. Check for correct polarity

5. Check for correct grounding

November 10th, 2022

ILT Panel Final Checks #1

PDI is squashed in SC control, needs more¹
CD62L more¹ TCCR4

SG bright vs average

PMA + CD69 are not in the dying CD3 pop

Might be better to do a spill in Zombie WFR control given gating

100

100

100

100

100

100

100

100

100

100

100

100

100

100

100

100