JMC: SEB/ Cord February 8th, 2023 PMA: Kitchen Sink 1:100 DNASET Specimen Status Location Conc Date **Tasks** Volume Ly Ly+Mon Total 3000 -34E46 15 EG NDOOG Adult 1/20/2023 9.26 m 5.63 2 6.87 13.74 (60/17 Inf-000374-4 HEU-10 -5 vials 868 6.68ml 2 4.34 598 11.96 NDOOG 10775 2.5 4.31 20.38pl 5.08 12.7 930c/p 4.225 Inf000374 2.5 3.69 4.84 11.72 12.1 9:04 an Than Indian m 9:39 3 stain 9:49 am Count stort (= 10% B cells 10:13 am Cells into incubater (6 nell plate 1 m/ml 5-6 mls) [5:18 pm] spining dam (6:20 madia manning) Cord = 9.225 million 15:43 pm stain ZMillionells 000 1.3, 00000 J: STPM Cells SEB (COUTA for CMs) 78:57 pm golgi bledl Adult 10.775 million 6 million sols Meder 1.5 % 3.69 = THOT HI + 593 HI RIO = 1.5 mm SEB 2.0 % 3.69 = 542/1 + 458/1 RIO = 20 El/L 10 /9/mL SEB (/pg/pm) = 10 pL/mL TBORF + 1 Monersin + 19 pt R10 > gots Mix -> 10pl of the mix into each tube. 2.225 Lym and > 2 vicils for undavid MI @ 1 will be SEB > 2 online 1 NM No 1 will be for Media > constant No 1 will be for Media > constant No 2008 2% 4.31 = 464 pl +536pl ?10 = 2.0\f6/mc + 10pl SEB + 10pl golgi max

wells SEB

Co 50 phrash cord in Adul In 000 0 a doses

10 doses 0 my 14/

(9:44 pm/ pulled from marbaker CD107A - 100 9:59 samples spin 1/6 Ime 10:06 45/c 08 cord -> 70 each scenstaineds 7 50 eads 10:32 4°C with nothing Single colors 4°C out single colors 37°c samples @ 10:45-711:15am 11:18 37°C Sample spin 1400 pm 6mon 11:29 4°C Sample stain 30 min (-> 11:59) _11:36 am scls start 1:10 11:42 Zambie single color cord 15 mm (> 11:57) 11:49 single colors 39°C(hot) stains 30 min (>12:19) 11:59 Zombie single color cord centrifug 1300 Hpm 8 mm (>12:07) 30- 72.5 - 37.5 Sels for 4° @ 12:00 pm -> [12:30 pm] 12:01 RBC lyse samples (> 12:04) LPC Samples spin (RBC lyxd) contrifuge samples 14 pm 6 m (=12:14) 12:08 -> 27 -> 37 -> 47 pm 12:17 12:21 PUBC type Single color hot 37° (> (>12:24)
12:27 REC lyse single color hot 37° (centridage 14 pm 6 m = (12:33) 12:30 \$ Single colors Gold 4℃ RBC lyse > (12:33) 13:31 3 . (1 11 (12:34) 37°C Stain 4°C 12:36 1:00 ,1:10 ,1:20 12:50 hot of cold single colors Ruclear Fix Perm

Samples wash 1 15 your 6 m, 12:59 second wash

Intracellular @ 1:09 > 1:49pm	
Sc's first Wash @ 1:27pm	
1:59 Sc/s away	
2:01 Samples Finish Wash 1	
Wash a not performed	
2:02 PFA JOHR added -> 4°C	
2:06 spikes single colors RT 40 mon -> (a:46)	
Done @ 9:04pm Aurora @ 9:40 am	
[QC Q. 9:42 am] event cate 140	
Passed	dims us bashl
// start @ 9.53am// from 250,000	dims of view pught
1/ Start @ 9.53am/ from 250000 Threshold to 150,000 eds FSC (other feture option is 1 FSC-A for	nuclear stains)
(other teters option is I FSC-A for	
Howard HOP ~ 1:30 see Egrandegto contour in	cord?
man CCRU - 40 exact	
2 sporadie see drops real, one on CD45RA cord	
Icos (stynal) Full? eleaghea? se CO3 =	t intra grile
blip co40L/ 1777	
	on the tolerial series
bo	Contraction series

		And UNSTAINED CONTROLS !!!		29 R8	28 R7		26 R4							19 84		83		1		-	12 V10		11 V7	10 v5	V3	9 43	1	7 UV16	6 U14	5 UV11	4 UV10	3 UV9	2 UV7	1 UV2
		III III		7						2002		9																						UV2
				APC/Fire 810	APC-Ffluor780	APC-R700	MIEXA FILIOF 64/	APC	PE-VI0770	PerCP-eF710	PE-Cy5	PE-CF594	PE	Spark blue 550	FITC	BV/86	BV750	BV711	BV650	BV605	BV570	BV510			Pacific Blue	Pacific Blue	BV421	BUV805	RHV737	CTOROL	RHV615	BUWEC2	AF	BUV395
			0000	CD38	מלון	CD107a	ILZ	CD39	TNFa	CD40L	CD25	CD26	PD-1	CD3	FoxP3	CCR6	IFNg	ICOS	CCR7	CXCR5		CD45RA			CD19	CD14	CD127	CACAS	CVCDO	CCK4	CD69	CD8	AF-UV6	CD62L
rippette draw volume /sample		R10 Media	Antibody Total	(CTIU)		(H4A3)	(MQ1-	(A1)			(M-A251)	(M-A261)	(EH12.2H7	(SK7)		(11A9)	(B27)		(G043H7)	(J252D4)		(Hi100)		(money)	(HIR10)	(MEET)	(ACTORE)	(JCb/CXCR3)		(161)	(FNS0)	(RPA-T8)		(DREG-56)
19.5		0.0		-		6.0	_				1			_	-							-		+	-	\mid	+	-	-	-		-	-	
						12																											1	
/sample	Brilliant Stain	Antibody Total			<1:2500>									1																				
62.7	50	15.7	1.6			1	1.4	3		1.2	1.2		1		1.5		+	1	2.5	٥	+					1,5	1.5	1		2	+	1	+	(0)
	100		3.2				4.1			2.4					u			2									5 3		-	2.5 5				@37C
61.4	50	14.4							1.2			1.2	1.2				1.2				0.7			2.0	2.0			X.		0.0	0./	100	1.2	SUMIN @4C
	100	28.8							2.4			2.4	2.4				2.4	0	200000		1.4			4	4		1	0		-	1.4		2.4	6
																						0.000						100	7				-	15 min
				-	_		_					200																						Fix/Perm
		13.2	1.5	1		1.2				0.5		0.5	0.5	5		1.5		1								0.5	2			1				@RT
	100	26.4	u	'		2.4		# i		1		_	1	10		w		2								1				2				2

First PermWash: Second Perm Wash:

2 ml NuclearWash 1500 rpm 6 min 2 ml NuclearWash 1500 rpm 6 min

Resuspend in 100 ul 0.4% PFA-PBS

tacos in cus inhancher

Simplified Protocol

Resuspend at 1E6 Lymint of R10, rest cells for <6> hrs in a 6-well plate Collect, count, aliquot cells 2.0E+6 Cells/ 5ml polystyrene tube Add antigers, <stim antibodies>, and CD107a in 1 ml R10 Cap and incubate at 37°C for 2 two hours
Add golgi block, and continue to incubate for 12 hours Thaw cells, DNAse, count.

800 ul of LiveDead mix (1:2500) @RT for 15min Wash 2 ml 5% PBS-FBS, spin 1300 rpm, 8min

Wash with 2 ml PBS, spin down 1300rpm 8min


Add to it HotStain mix, Incubate @37C for 30 min Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Add ColdStain mix, incubate @ 4C for 30min Add 300-500 ul 1x RBC Lysis for 3 minutes Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

1 ml Nuclear FixPerm, incubate @ 4C for 30min

(vortex every 10 minutes)
2 ml NuclearWash 1500 rpm 6 min
2 ml NuclearWash 1500 rpm 6 min

Add Intracellular Stain, incubate @ RT for 40min
First PermWash: 2 ml Nuclear
Second Perm Wash: 2 ml Nuclear



		And UNSTA			29	28	27	26	25		-		21	20	19	18		16	15	14	ti	12	=	:	10		9	∞	7	6	5	4	w	2	.	==
		And UNSTAINED CONTROLS !!!		172	88	87	R6	R4	R2	RI	813	810	888	86	84		Ra	VIS	VIA	V13	VII	V10	*		V5	V3	V3	V1	UV16	U14	UV11	UV10	6AN	UV7	240	Filter Single color (ul) Ref ctrl type name
				APC/Fire 810	APC/Fire 750	Lombie NIR	APC-R700	Alexa Fluor 64/	APC	PE-FIFEXIU	PerCP-eFlour710	PE-Cy5	PE-CF594	PE	Spark blue 550	FITC	BV786	BV750	BV711	BV650	BV605	BV570	BV510	BV480		Pacific Blue	Pacific Blue	BUACO	BUVISI	Tagand	BUV615	BUV563	BUV496	AF		Fluorochrome
				CD38	4-1BB	L/D	CD107a		CD39	OX40	CD40L	CD25	CD26	PD-1	CD3	FOXP3	CCR6	IFNg	icos	CCR7	CXCR5		CD45RA		10000	CD19	(017/	CD4	CXCR3	Adz	CCR4	CD69	CD8	AF-UV6	CD62L	Marker
Pippette draw volume /sample		R10 Media	Antibody Total	(НІТ2)			(H4A3)		(A1)			(M-A251)	(M-A261)	(EH12.2H7	(SK7)		(11A9)	(B27)		(G043H7)	(J252D4)		(Hi100)		(Crain)	(MDEZ)	(A019D5)	(SK3)	(1C6/CXCR3)	(80)	(161)	(FN50)	(RPA-T8)		(DREG-56)	Clone Vial or #
19.5	14.5	4 00	5				6.0											1	1							-			-	-		-	+	+		stimili
	-						12																													
Pippette draw volume /sample	Brilliant Stain	Antibody Total			-T-2000	<1.2500>					1																								Total State of the last of the	15 min
ne 62.7	50	15.7	1.6			1	1	1.4		1	7.7	7.1		1	-							-	i	2					-			+	+	-	1000	min
7	100	7 31.4	3.2					4.7		-	2 2.4					1.5	-		1 2	2.5							1.5	2		1.1					(Corc	SUITIO
64.5	50	17.5		1.2					1.2	1.2			1.2	1.2				1.2				0.7			2.0			8	7.0		0.5	0.7		1.2		2 30min @4C
	100	35.0		2.4					2.4	2.4			2.4	2.4				2.4				1.4				4		3	1		1	1.4		2 2.4		@4C 2
																						2016	<0.8>	Stren				2.0	7.			**		43	Fix/Perm	16 min
58.7	50	11.7					1.2				0.5		0.5	0.5	υī		1.5		1								0.5				1				COLUMN TO	Stain 40 min
1	100	23.4	I				2.4				1		1	1	10	68	w	1	2				1316	-	+	-	100			1000		681	-		100	

> (6)

2.69.5

Simplified Protocol

Thaw cells, DWAse, count.

Requisipend at 1EE Lymmin of R10, rest cells for <6> hrs in a 6-well plate Collect, count, aliquot cells 2.0e+6 Canla fin polystyrene lube. Add antigens, salim antibodiesy, and CO107a in 1 ml R10
Cap and incubate at 31°C for 2 two hours
Add poligi block, and continue to incubate for 12 hours

Add to it HotStain mix, incubate @37C for 30 min Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Wash with 2 ml PBS, spin down 1300 pm 8min 800 ul of LiveDead mlx (1:2500) @RT for 15min Wash 2 ml 5% PBS-FBS, spin 1300 rpm, 8min


Add ColdStain mix, incubate @ 4C for 30min Add 300-500 ul 1x RBC Lysis for 3 minutes Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

1 ml Nuclear FixPerm, incubate @ 4C for 30min

First PermWash: Second Perm Wash: (vortex every 10 minutes)
2 ml NuclearWash 1500 rpm 6 min
2 ml NuclearWash 1500 rpm 6 min

Add Intracellular Slain, incubate @ RT for 40min
First PermWash:
2 ml NuclearWash 1500 rpm 6 min
Second Perm Wash:
2 ml NuclearWash 1500 rpm 6 min

Resuspend in 100 ul 0.4% PFA.PBS
Cap tubes, wrap rack in foil, store at 4*C



3 MM O O O O Cgf

0 0 0 0 Cgf

3.2 (2 md/s PMH) OSE 8 Adil 000 L

1.00 EN 155

And UNSTAINED CONTROLS !!!			R8 APC/Fire 810	R7 APC/Fire 750	R6 Zombie NIR		R2 AlexaFluorb4/	R1 Alexa Bloom 647			g.		99			Ale					BV605	V7							BUV737	DATA BOATA			В	UV2 AF	SEEVUB SEEVUB	Hiter Single color (ar) name
		1	CD38	CD27	U/D	CDTO/a	CD407s	Va7.2	CD16	PD1			CD26	NKG2D	CD3	hCD1d	CCR6	IFNY	CD7	CCR7	CD56		CD45RA					CD4	CXCR3	V82	CCR4	CD69 (Ĭ	CD62L (DR	
Pippette draw volume /sample	R10 Media	Antibody Total	(НІТ2)	(0323)		(DAMO)	(EAAA)				(MAB11)	(M-A251)	(M-A261)		(SK7)		(11A9)	(B27)					(Hi100)	(REA631)	(нів19)	(MSE2)	(A019D5)	(SK3)	(1C6/ CXCR3)	(B6)	(161)	(FN50)	(RPA-T8)		(DREG-56)	
19.5	14.5	6.0		1		410	6.0																													
	116	6	5			_	48																	-						-						
Pippette draw volume /sample	Brilliant Stain		Antihody Total			<1:2500>		<: >()								()<:>()																				
le 66	y y	5	19.0	1.6	2							7.1	17	1 3	4		1.0	2 7	+	-> -	1.0	10		7	2		T.O	1.5	1		1.0	20				The second second
	400			12.8	16			100000				0.1	9.6	9.6	13		-	10		00 0	000	o o		TO	16		7.7	12	13		1	16				
60.7			13.7	8	5,			4.4	1 2	0.7	1.5					1.2	c > (1.2)							0.7	1:0	20	20		1:0	12	0.7	0.0	0.5	0.7	1.2	
	100	400	109.6					0.0	9.6	5.6	12.0					9.6								5.6		16.0	16.0			10.4	5.6		4.0	5.6	9.0	2
			6										350																							
26.5	n n	50	9.5						0.5		0.5	2.5	0.5		0.5	0.5			1.5		1								0.5		0.5		1			
		400	76						4.0		4.0	20.0	4.0		4.0	4.0			12.0	TO SERVICE SER	8.0				5.0				4.0		4.0		8.0			SIGNATURE

Simplified Protocol

Thaw cells, DNAse, count.

Collect, count, aliquot cells 3.0E+6 Cells R10 / 5ml polystyrene tube
Bring volume upto 1 ml R10, add 2 ul PNAXCtrl and CD107a
Cap and incubate at 37°C for 6 hours

Wash with 2 ml PBS, spin down 1300rpm 8min 800 ui of Livebead mix (12500) @RT for 15min Wash 2 ml 5% PBS-FBS, spin 1300 rpm, 8min Add Holstain mix, incubate @37C for 30 min Wash 2 ml 5% PBS-FBS 1400 rpm, 6 min

h 2 ml 5% PBS-FBS 1400 fpm, o min

Add Tetramers, incubate @RT for 10 min Add CodeStain mix, incubate @ 4C for 30min Add 300,500 ul 1x RBC Lysis for 3 minutes Wash 2 ml 5% PBS-FBS 1400 rpm, 6min 300 ul BD FixPerm, incubate @ 4C for 20min

(vortex every 10 minutes)

Wash: 1 ml PermWash 1500 rpm 6 min erm Wash: 1 ml PermWash 1500 rpm 6 min

Add Intracellular Stain, incubate @ RT for 40min First PermWash: 2 ml PermWash 1500 rpm 6 min

Resuspend in 70 ul 0.4% PFA-PBS Cap tubes, wrap rack in foil, store at 4°C

513.6 8 GIV in soul

PMA: Kitchen Sink Initial Unmixo (splitting aggured into 2 for unmixing) TNFa on Pe-Cy7 nearly off scale is (for adult) No dear signal for IL-17a (vide occertants bands

AC040L

wait... intracellular spille

90×109×10 90×105

