

August 67, 2021

HEU: P2 + P4b $\gamma\delta$ Cord Optimization #2

Specimen	Status	Location	Conc	Date	Tasks	Volume	Ly	Ly+M	Total	1.5E+6	2E+6
3004 <i>NY054</i>	<i>Adult</i>		<i>1047</i>	<i>11/12/15</i>	<i>Pub tests</i>	<i>3 ml</i>	<i>3.56E4</i>	<i>4.68E4</i>	<i>10.68</i>	<i>420pl</i>	<i>941 sc's</i>
3804 <i>CQ1047</i>	<i>Cord</i>					<i>3 ml</i>	<i>1.30E4</i>	<i>1.51E4</i>	<i>4.53</i>		<i>1.54mb</i>

June	Ex vivo Gd-p2	Surface μL x	Intra μL x
BV421	VD2	1	
BV510	CD3	1.5	
BV650	NKG2D	2.0	
Alexa 488 (FITC)	dg9	-	3
PerCPeF710	CD56	1	
PE	PD1	1.5	
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	VD1	1.0	
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$	PBS	<i>13.4</i>	17.5

June	Ex vivo Gd-p2	Surface μL x	Intra μL x
BV421	VD2	1	
BV510	CD3	1.5	
BV650	NKG2D	<i>2.2</i>	
	Streptavidin	<i>1.5</i>	
Alexa 488 (FITC)	dg9	-	3
PerCPeF710	CD56	1.0	
PE	PD1	1.5	
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	VD1	1.0	
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$	PBS	<i>11.9</i>	17.5

Thaw start @ 10:08 am
 Spin @ 10:18 60pl DUA (not very many cord...)
 so stain @ 10:37

*Guava being hellu uncooperative this Am)
 very few monocytes in this cord blood sample.

Plating tubes @ 11:10 am

★ Made new scs. to be safe ★

L/D @ 11:42

Scs @ 11:54 \rightarrow 12:19 \rightarrow Done @ 12:35 pm

Surface @ 12:19

RBC lysis 12:42 \rightarrow 12:46

Biotin @ 12:59 pm +15 \rightarrow 1:14 ✓

Fix Perm @ 1:02 \rightarrow 1:12 \rightarrow 1:22

strip Fix Perm @ 1:26 pm \rightarrow 1:36 \rightarrow 1:46

Intra cell lysis @ 1:43 pm \rightarrow 2:23 pm

@ 2:03 pm \rightarrow 2:43 pm

Done @ 14:56 pm

3pl - 3ml Horizon

4x8 = 3.2 3.5ml - 7pl Aqua

Streptavidin \rightarrow before or
 2ml + 1.5pl

after RBC lysis
 "for this round"

Strip 14ml inc @ RT
 wash w/ 2ml 5% FBS PBS
 (Δ vs 1 ml or PBS?)

August , 2021

Testing on an adultHEU: P2 + P4b $\gamma\delta$ Cord Optimization #2

Hao-Ting & November	Ex vivo Gd-p4b	Surface μL x	Intra μL x
BV421	PD1 ✓	2	
BV510	Aqua L/D	1:500	
BV650			
Alexa 488 (FITC)	Perforin (dG9)	-	- 2
PerCPeF710	CD3	1.5	
PE	GZMb	-	- 1.5
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	CD56	1.3	
APC Fire750	VD2	1.2	
20.5 $\mu\text{L}/\text{rxn}$;	PBS	12	17.0

Switchup	Ex vivo Gd-p4b	Surface μL x	Intra μL x
BV421	PD1	2	
BV510	Aqua L/D	1:500	
BV650			
Alexa 488 (FITC)	Perforin (dG9)	-	- 2
PerCPeF710	CD3	1.5	
PE	GZMb	-	- 1.5
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	VD2	1.3	
APC Fire750	CD56	1.2	
20.5 $\mu\text{L}/\text{rxn}$;	PBS	12	17.0

(13.9) - November

Current Version	Ex vivo Gd-p4b	Surface μL x	Intra μL x
BV421	PD1 ✓	2	
BV510			
BV650	CD56	1.8	
Alexa 488 (FITC)	Perforin (dG9) ✓	-	- 3 16.5
PerCPeF710	CD3 ✓	1.5	
PE	GZMb ✓	-	- 1.5 8.25
PE Dazzle			
PE Vio770	NKG2A ✓	0.6	
APC	V82	1	
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$;	PBS	13.6	16 88

13.1

Test Aqua	Ex vivo Gd-p4b	Surface μL x	Intra μL x
BV421	PD1	2	
BV510	Aqua L/D	1:500	
BV650	CD56	1.8	
Alexa 488 (FITC)	Perforin (dG9)	-	- 3
PerCPeF710	CD3	1.5	
PE	GZMb	-	- 1.5
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	V82	1.0	
APC Fire750			
20.5 $\mu\text{L}/\text{rxn}$;	PBS	13.6	16

Test Aqua + CD16	Ex vivo Gd-p4b	Surface μL x	Intra μL x
BV421	PD1	2	
BV510	Aqua L/D	1:500	
BV650	CD56	1.8	
Alexa 488 (FITC)	Perforin (dG9)	-	- 3
PerCPeF710	CD3	1.5	
PE	GZMb	-	- 1.5
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	V82	1	
APC Fire750	CD16	1	
20.5 $\mu\text{L}/\text{rxn}$;	PBS	12.6	16

68

7

0

0 1

August 6-7, 2021

HEU: Innate-like T cell Panels
Tetramer Panels

	Ex vivo MAIT-p1a	Surface μL x	
V450 BV421	Va7.2	2	
V525 BV510	Aqua L/D	1:500	
V670 BV650			
B530 FITC	CD4	(1.5)	
B710	CD45RA	2	
PerCPeF710	CD161	0.5	
Y590 PE			
Y615 PE Dazzle			
Y780 PE	CD3	0.5	
Vio770			
R670 APC	CCR6	(2.5)	
R780	CD8	(1.2)	
APCFire750			
20.5 $\mu\text{L}/\text{rxn}$;	PBS	10.3	

Concurrent Staining (20 minutes, 4°C)
(1:1000)

Separate Staining (20 minutes tetramer RT, 20 minutes 4°C)
(1:1000)

Tetramer Test A1	Ex vivo MAIT-tet	Surface μL x	
V450 BV421	CD4	1.5	
V525 BV510	Aqua L/D	1:500	
V670 BV650			
B530 Alexa488 (FITC)	6FP/ 5OP-RU		
B710	CD45RA	2	
PerCPeF710	CD161	0.5	
Y590 PE			
Y615 PE Dazzle			
Y780 PE	CD3	0.5	
Vio770			
R670 APC	(CCR6)		
R780	CD8	1.2	
APCFire750			
20.5 $\mu\text{L}/\text{rxn}$;	PBS		

Tetramer Test A2	Ex vivo MAIT-tet	Surface μL x	
V450 BV421	CD4	1.5	
V525 BV510	Aqua L/D	1:500	
V670 BV650			
B530 Alexa488 (FITC)	6FP/ 5OP-RU		
B710	CD45RA	2	
PerCPeF710	CD161	0.5	
Y590 PE			
Y615 PE Dazzle			
Y780 PE	CD3	0.5	
Vio770			
R670 APC	(CCR6)		
R780	CD8	1.2	
APCFire750			
20.5 $\mu\text{L}/\text{rxn}$;	PBS		

Table 1: Summary of Data

Category	Value 1	Value 2	Value 3
A	10	20	30
B	40	50	60
C	70	80	90
D	100	110	120
E	130	140	150
F	160	170	180
G	190	200	210
H	220	230	240
I	250	260	270
J	280	290	300

Table 2: Detailed Data

Category	Value 1	Value 2	Value 3	Value 4
A	10	20	30	40
B	40	50	60	70
C	70	80	90	100
D	100	110	120	130
E	130	140	150	160
F	160	170	180	190
G	190	200	210	220
H	220	230	240	250
I	250	260	270	280
J	280	290	300	310

August , 2021

HEU: P2 + P4b $\gamma\delta$ Cord Optimization #2

Flow panels set up @ 10:33 am ✓

- adjusted gate to cut out slumpy curve on all...



Run start @ 10:48 am ✓ ↑ 580 FSC

Biotin'd NRG2D does look better.... pd1 balls squashed?
CD3 -

★ Clot ★ 11:06 am ✓

clot resolved @ 11:11?

★ water tube had too cells ★

② @ stop acquire, how long tube to fall back in?

- Does this happen?

B running dirty 3500 cts/sec

★ BD heavy non specific ★

Done @ 11:50 am

Culture Set up:
1.5 EC LTM/ml

D₀
D₁₀
D₁₂