

November 9-10<sup>th</sup>, 2021

HEU: Ex-vivo  $\gamma\delta$  CBMC panels #10 & ILT novo 3

$$\text{Inf } 121 \quad 21 = 356 + 344$$

$$\text{Inf } 124 \quad 21 = 500 + 200$$

$$\text{Inf } 171 \quad 22 = 295 + 1905$$

$$\text{Inf } 129 \quad 22 = 472 + 228$$

$$\text{NY062 } 434 + 246$$

$$\begin{array}{r} \nearrow 20 \\ 10 \text{ PMA} \\ \times 3 \\ \hline 3 \text{ ml} \\ +16\% \quad 31.5 \end{array} \quad \begin{array}{l} a \text{ ctrl} \\ 27 \text{ ml} \\ 224 \text{ ml} \end{array}$$

Used old protein transport inhibitor  
(2019) as not enough

7:21 pm  $\rightarrow$  spindown

Tet  $\frac{1}{10}$  @ 7:36

everything vortexed @ 8:08 pm

Tetramers @ 8:35 pm

Adult PBMC libraries setup

NY062  $\rightarrow$  stimulus PMA  
 $\rightarrow$  ctrl w/o PMA

Te

Mass single color spin @ 9:56  
37°C @  $\rightarrow$  10:04 pm

100  $\mu$ l Fix Perm

4°C @ 10:25 pm  $\rightarrow$  10:55

Mass Fixperm @ vortex @ 22:37:47

Rbc lysis @ 11:02 pm

Tet Fix Perm @ 11:26 pm

Intra @ 00:13 pm

Sc's intra @ 23:50

black - PMA?  
blue - ctrl

4,100,000  $\div$  30  $\rightarrow$  1.3K cells/tube  
in ~100  $\mu$ l volume  
 $\rightarrow 1.3 \times 10^6$

Not confident anything is in there, but  
we shall see tomorrow if  
they are, or we just wasted  
time, effort & material



#	Filter	Single color (u)	Ref crt	Fluorochrome	Marker	During stimulation	Tetramer (RT) for 40 min	37°C (RT) for 30min	Surface staining 40C for 30min	Lysing, then CytoFix/Perm	RT for 40min Universal Intra mix (10x diluted)	
1	UV2	1.25	Cells	BUV395	CD62L (DREG-56)							
2	UV7	1.25	Cells	BUV496	CD38 (RPA-18)				1.0	6.0	0.25	1.5
3	UV9	1.25	Cells	BUV563	CD69 (REG-50)				1.0	6.0	0.25	1.5
4	UV10	1.25	Cells	BUV615	CD4 (IG1)			2.0	12			
5	UV11	1.25	Cells	BUV661	Vα2 (B6)				1.0	6.0		
6	UV14	1.25	Cells	BUV737	CD38 (156/CD38)			2.5	15			
7	UV16	1.25	Cells	BUV805	CD4 (S3)				1.5	9.0	0.25	1.5
8	V1	1.25	Cells	BUV421	CD127 (A01D05)	0.8			1.0	6.0		
9	V3	1.25	Cells	Pacific Blue	CD14 (M5F2)				2.0	12.0		
10	V5	1.25	Cells	Pacific Blue	CD19 (HB19)	0.8		1	6		0.5	3
11	V7	1.25	Cells	BUV480	CD161 (REA631)				1.0	6.0		
V10			Cells	BUV510 (dim)	CD45RA (HI100)							
V11			Cells	BUV605								
V13			Cells	BUV711								
12	V14	1.25	Cells	BUV750	IFNγ (B27)						1.5	9
13	V15	1.25	Cells	BUV786	CCR6 (11A9)			1.8	10.8			
14	B2	1.25	Cells	Alexa488	hCD1d		<1.750>					
15	B3	1.25	Cells	Spark blue 550	CD3 (SK7)				1.2	7.2	0.25	1.5
B4			Cells	PE								
16	B6	1.25	Cells	PE-CF594	CD26 (M-A261)			1.2	7.2			
17	B8	1.25	Cells	PE-CF5	CD25 (M-A251)				1.2	7.2	0.25	1.5
18	B10	1.5 (<10.4)	Cells	PerCP-Cy5.5	TMβa (MAB11)						2.5	15
B13			Cells	PE-wo770								
19	R1	1.25	Cells	APC	IL-4						1.5	9
20	R2	1.25	Cells	Alexa Fluor 647	hMFI SOP-RU		<1.500>					
21	R4	1.25	Cells	APC-R700	CD107a (H4A3)	6.0						
22	R6	1.2500	Cell	Zombie NIR	Viability							
23	R7	1.3	Cells	APC/Fire 750	CD27 (O323)			2	12			
24	R8	1.3	Cells	APC/Fire 810	CD38 (H12)			1.6	9.6			
Total volume						7.6		12.1	72.6	14.1	84.6	7.25
unstimulated Cells								50	50	50		44
24 color CBMC								57.1	59.1	52.3		

Protocol		
Aliquot cells		
PMA-Ionomycin stimulation		
Golgi blockers		
PBS 2ml	1400rpm, 6min	
Live/Dead @ RT 15min		
PBS+FB5 2ml	1400rpm, 6min	
Tetramers @ RT 40min		
PBS+FB5 2ml	1400rpm, 6min	
Surface	CCR4 CD26	
37°C 30min	CCR3 CD27	
	CD161 CD38	
PBS-FBS 2ml	CCR6	
Surface mix @ 40C 30min		1400rpm, 6min
Lysing solution 300-500ul		
PBS-FBS 2ml		1400rpm, 6min
Cytofix/perm @ 40C 20min		
Wash twice PermWash 1 ml		1500rpm, 6min
Intra Staining @ RT for 40min		
Wash once w/ 2ml		1500rpm, 6min
Resuspend 100 ul PFA-PBS		

CCR3 in MFI  
+ new CD4  
+ new CD4

RNA & cells got subcloned?  
Keep record  
on 12/11/21  
good thing we  
did

20:50  
11:45

Hot @ 8:56 → 12/11 (10min HT)  
Cold @ 9:15 → 12/11 (10min HT)  
25 → 34 → 16:04 PM  
Unstained (check @ 9:40)  
mass spin @ 9:50

Did I add them? (PDI, CD56, M1625?)





*PBMC control*

#	Filter	Single color (ul)	Ref ctrl	Fluorochrome	Marker	During stimulation	Tetramer (RT) for 40 min	37°C (RT) for 30 min	Surface staining 40°C for 30 min	Lysing, then Cytofix/Perm	RT for 40 min Universal Intra mix (10x diluted)	
1	UV2	1.25	Cells	BUV395	CD62L (DMEG-56)							6
2	UV7	1.25	Cells	BUV496	CD8 (RPA-78)				1.2	7.2	0.25	1.5
3	UV9	1.25	Cells	BUV563	CD69 (RM50)				1.0	6.0	0.25	1.5
4	UV10	1.25	Cells	BUV615	CCR4 (IG1)			2.0	1.0	6.0	0.25	1.5
5	UV11	1.25	Cells	BUV661	Vα2 (B6)				1.0	6.0		
6	UV14	1.25	Cells	BUV737	CXCR3 (156/CXCR3)			2.5	1.5	9.0	0.25	1.5
7	UV16	1.25	Cells	BUV805	CD4 (SK3)				1.5	6.0		
8	V1	1.25	Cells	BUV421	CD127 (A019D5)	0.8			1.0	6.0		
9	V3	1.25	Cells	Pacific Blue	CD14 (M5E2)				2.0	12.0		
10	V5	1.25	Cells	Pacific Blue	CD19 (HB19)			1	2.0	12.0	0.5	3
11	V7	1.25	Cells	BUV480	CD161 (REA631)	0.8			1.0	6.0		
	V10	1.25	Cells	BUV510 (dim)	CD45RA (HI100)							
	V11	1.25	Cells	BUV605								
	V13	1.25	Cells	BUV711								
12	V14	1.25	Cells	BUV750	IFNγ (B27)						1.5	9
13	V15	1.25	Cells	BUV786	CCR6 (11A9)			1.8	10.8			
14	B2	1.25	Cells	FTIC	Vα24a18				0.5	3.0		
15	B3	1.25	Cells	Spark blue 550	CD3 (SK7)				1.2	7.2	0.25	1.5
16	B6	1.25	Cells	PE								
17	B8	1.25	Cells	PE-CF594	CD26 (M-A261)			1.2	7.2			
18	B10	1.25	Cells	PE-Cy5	CD25 (MAA251)				1.2	7.2	0.25	1.5
	B13	1.25	Cells	PerCP-Cy5.5	TFNα (MAB11)						2.5	15
19	R1	1.25	Cells	APC	IL-4							
20	R2	1.25	Cells	Alexa Fluor 647	hMA1 6-FP						1.5	9
21	R4	1.25	Cells	APC-R700	CD107a (HA43)	6.0	<1:500>					
22	R6	1.2500	Cells	Zombie NIR	Viability			2	12			
23	R7	1.3	Cells	APC/Fire 750	CD27 (O323)			1.6	9.6			
24	R8	1.3	Cells	APC/Fire 810	CD38 (HR2)							
Total volume						7.6	12.1	72.6	14.6	87.6	7.25	44
unstimulated cells												
24 color CHMC							Brilliant Stain: 50	Draw Volume: 57.1	50	59.6	50	52.3

*Rumins PMA ITI wa extra first*

*on the control, right well group (F) unstimulated control*

*6x5=30*  
 $\frac{4-36}{4} = 6 \neq 7$   
 $\frac{4-36}{4} = 6 \neq 7$   
 $\frac{4-36}{4} = 6 \neq 7$

Protocol	
Aliquot cells	
PMA-ionomycin stimulation	
Golgi blockers	1400rpm, 6min
PBS 2ml	
Live/Dead @ RT 15min	
PBS-FBS 2ml	1400rpm, 6min
Tetramers @ RT 40min	
PBS-FBS 2ml	1400rpm, 6min
Surface	CCR4 CD26
37°C 30min	CCR3 CD27
	CD161 CD38
PBS-FBS 2ml	1400rpm, 6min
Surface mix @ 40°C 30min	
Lysing solution 300-500ul	
PBS-FBS 2ml	1400rpm, 6min
Cytofix/perm @ 40°C 20min	
wash twice PermWash 1 ml	1500rpm, 6min
Intra Staining @ RT for 40min	
wash once w/ 2ml	1500rpm, 6min
Resuspend 100 ul PFA-PBS	

*Stim +/-*

*Stim +/-*

*Stim +/-*



No CD56, PD1 or NKG2D got added last night.

CCR3 ≠

my Va24jal8 stain at 0.5 not present.

hmr1 + IL-4 really ugly q. w. h

hmr1 vs uclz.

IFN $\gamma$  / TNF $\alpha$  not here?

did I use TNF $\alpha$ ?

Yeah ~ IL-4 ~ Va24jal8

DTR- PBMC -  $\frac{2}{b}$  - 2021

11/11/21

Questions: Blue Badl  
PMA ~~ctrl~~ ? (did it even PMA?)

- are there any cells there?

Synergy ≠ robust

(no apparent PMA well issue, no stim)  
Vortex everytime w/dilution ~ 200 events  
there are cells there ~ 220 events.

CD62L → run at, included "dying"  $10^3 \frac{1}{2} \rightarrow 10^6$  ... not trust manually

CD8 → looks PMA-like

CD69 looks activated ✓  
+ no CD62L → suggest you did work!

You can remove tubes! ✓

Blackl → PMA

Quite a few odd events  $10^6$ , gains not set?

- badl CD19 looks odd -

important pop markers do 50,000 (CD161)

- yup these stims worth ed... wtf happened other wells?


brightest CD3's (not the downreg stim ones)

no overwork until total cell collected

NKG2D run to 50k (not that bright).

TNFA ≠ there? (is gas for IFN $\gamma$ ?)

- We want the highest PD1 expressors, exclude it to 50k -  $\rightarrow$  is there ✓



Does it make sense?  
CCR6 also in this context!

Ref control next (moment of truth)

work needed it!

Start on 2<sup>nd</sup> set @ 9:55pm

DTR-PBNCb - -2021

F  
chl

May not have enough vD2's,  $\frac{1}{2}$  blue segments across spectrum

Some specimens, shift to the left you gain significance values (cleared average spread?)

Full color →

Alex Mei

Distinguishing AF from highest expressors when gating.

highest CD3%  $\geq$  cd/y's?

"CD25 just the tip"

It's the brightest that show overcompensation/unmixing issues,  
the average are generally fine.

... only vD2? ...

forget to <sup>I or whatever</sup> gate adjust  
panel

11:53pm