

March 7th, 2023

1000 counts:

1:100

2023\_ILT\_01

2

R10

| Specimen        | Status               | Location | Conc                  | Date    | Tasks                    | Volume | Ly   | Ly+Mon | Total | 3E+6               | R10         |
|-----------------|----------------------|----------|-----------------------|---------|--------------------------|--------|------|--------|-------|--------------------|-------------|
| Info189<br>a-3  | ♂ HU<br>RTL8         |          | 11.6 $\mu$ l<br>(12)  |         | bloody<br>1              | 2      | 5.06 | 6.76   | 10.12 | 395 $\mu$ l<br>(5) | 205 $\mu$ l |
| Inf161-8<br>a-3 | ♀ HEV-6<br>SBKC      |          | 20.13 $\mu$ l<br>(18) |         | not bloody<br>3 - clumps | 2      | 5.01 | 8.09   | 10.02 | 349 $\mu$ l<br>(5) | 201 $\mu$ l |
| Inf180-9<br>a-3 | ♂ HEV-<br>hi<br>sqWB |          | 20.7 $\mu$ l<br>(15)  |         | not bloody<br>2          | 2      | 3.85 | 5.27   | 7.7   | 519 $\mu$ l<br>(4) | 81 $\mu$ l  |
| ND050           | Adult<br>Normal      |          | 21.2 $\mu$ l<br>15EG  | 1/13/23 | 1                        | 2      | 4.64 | 6.35   | 9.28  | 431 $\mu$ l<br>(5) | 169 $\mu$ l |
| ND006           | Adult<br>Sc's        |          | 23.0 $\mu$ l<br>15EG  | 1/20/23 | 2                        | 2      | 3.55 | 4.93   | 7.1   | 563 $\mu$ l<br>(3) | 47 $\mu$ l  |

9:31 am cord thaw

9:44 am adult DNAse

Count stain @ 10:13 AM

10:28 am count

Guava was having issues after being left on for 4 days  
and plate overfilled on guava...  
and pellets looked solid... Guava's count is off?

→ Run @ 2EG cells

(See what turns upon Aurora)  
~ 3 EG cells

11:41 am Incubation  
Start, lower rack  
@ 37°C.  
2 Million Lym/ml  
2 ml

→ 5:41 pm out

Ab spin @ 4:34 pm

@ 5:31 pm

Sample spin @ 5:57 pm

46 @ 6:18 pm

Sc's 6:34

caliquels → fudge

Hot samples @ 6:58 pm → 7:28 pm

Sc's spin @ 7:09 pm

Hot sc's @ 7:35 pm → 8:05 pm

Cold sc's @ 8:00 → 8:30 pm spin @ 8:34 pm

Samples to 3°C @ 8:00 pm Pause!!!

8:23 hot sc's paused pre-fix

8:29 tetramers

8:40 pm Va7.2 / Va24

(check background, 1.2 30 min RT vs 20 min RT?)

Sc's FixPerm w/ 200  $\mu$ l @ 8:51 → 9:01 → 9:11 pm

Samples cold @ 9:32 pm → 10:02 pm

Samples post RBC lyse spin @ 10:11 pm

10:20 pm intrac sc's → 11:00 pm

10:25 FixPerm → 35 → 45

2  $\mu$ l/ml

\*  
400  $\mu$ l R10 400  
x 7.2 x 15  
800 2000  
28000 4000


2880  $\mu$ l & 144  $\mu$ l  
FMT

6000  $\mu$ l & 30  $\mu$ l  
goig: black

4  
15  
x 8

12.0 mls

5  $\mu$ l in 12.5 mls



CERT → 14  $\mu$ l into intracellular,  
rest for sc.

1.2 20.5  
x 10 x 10  
12  $\mu$ l 205  $\mu$ l  
- 12  
193  $\mu$ l

going forward

1.2 for abs  
Bpl 1:100 let

Final sc spin & 2nd <sup>sample</sup> Perm Wash @ 11:07pm

Sc's done @ 11:18pm

Samples intake @ 11:27 (w/ appropriate Vol. 2/Vol 24 spikes)  
→ 12:07pm

Done @ 12:24am

2015  
01/27  
01/27

01/27

| #                            | Filter | Single color (ul) | Ref ctrl | Unstained ctrl | Fluorochrome | Matter | Clone | Vial Lot # | During  | 14 | U/D 15 min | Tetramer 40 | Host        | 14         | 14   | RBC Lys | Splend | 14 |
|------------------------------|--------|-------------------|----------|----------------|--------------|--------|-------|------------|---------|----|------------|-------------|-------------|------------|------|---------|--------|----|
| 1                            | UV2    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            | stim II |    | (RT)       | min @ RT    | 30min @ 37C | 30min @ 4C | 16.8 | then    | @ RT   |    |
| 2                            | UV7    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.2        | 16.8 |         |        |    |
| 3                            | UV9    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.7        | 9.8  |         |        |    |
| 4                            | UV10   | ND050 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.5        | 7.0  |         |        |    |
| 5                            | UV11   | ND050 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 2.0        | 28   |         |        |    |
| 6                            | UV14   | ND050 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.7        | 9.8  |         |        |    |
| 7                            | UV16   | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.3        | 18.2 |         |        |    |
| 8                            | UV1    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21   |         |        |    |
| 9                            | UV3    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21   |         |        |    |
| 10                           | UV5    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 2.0        | 28.0 |         |        |    |
| 11                           | UV7    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 2.0        | 28.0 |         |        |    |
| 12                           | UV10   | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.7        | 9.8  |         |        |    |
| 13                           | UV11   | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.0        | 14   |         |        |    |
| 14                           | UV13   | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.0        | 14   |         |        |    |
| 15                           | UV14   | ND050 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.0        | 14   |         |        |    |
| 16                           | UV15   | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.0        | 14   |         |        |    |
| 17                           | B2     | ND050 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21   |         |        |    |
| 18                           | B3     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21   |         |        |    |
| 19                           | B4     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.2        | 16.8 |         |        |    |
| 20                           | B6     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21   |         |        |    |
| 21                           | B8     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.2        | 16.8 |         |        |    |
| 22                           | B10    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.2        | 16.8 |         |        |    |
| 23                           | B13    | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21.0 |         |        |    |
| 24                           | R1     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.7        | 9.8  |         |        |    |
| 25                           | R2     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 1.5        | 21.0 |         |        |    |
| 26                           | R4     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.5        | 7.0  |         |        |    |
| 27                           | R6     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.5        | 7.0  |         |        |    |
| 28                           | R7     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.5        | 7.0  |         |        |    |
| 29                           | R8     | ND006 chl         |          |                | AF           | AF-UV6 | SK11  |            |         |    |            |             |             | 0.5        | 7.0  |         |        |    |
| And UNSTAINED CONTROLS III   |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| R10 Media                    |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| Pipette draw volume / sample |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| Antibody Total               |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 6.0 84                       |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 14.5 203                     |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| Pipette draw volume / sample |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 19.5                         |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| Antibody Total               |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 19.0 266                     |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 50 700                       |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 59.5                         |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 175.0                        |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 700                          |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 9                            |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 126                          |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 50                           |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 700                          |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |
| 56                           |        |                   |          |                |              |        |       |            |         |    |            |             |             |            |      |         |        |    |

ND006 activated  
→ RBC

• No CD3 or PDI in hka spile  
• No surface ccr27 (1/2 dose added only in hka)

20.5  
-24  
181  
Simplified Protocol

Thaw cells, DNase, count.  
Collect, count, aliquot cells 2.0E6 cells R10 5ml polypropylene tube  
Bring volume up to 1 ml R10, add 2 ul PI/ACGI and CD107a  
Cap and incubate at 37C for 8 hours

Wash with 2 ml PBS, spin down, 1300rpm 8min  
800 ul of Live/Dead mix (1:2500) @RT for 15min  
Wash 2 ml 5% PBS-FBS, spin 1300 rpm, 8min  
Add HostStain mix, incubate @37C for 30 min  
Wash 2 ml 5% PBS-FBS 1400 rpm, 6 min  
Add Tetramers, incubate @RT for 10 min  
Wash 2 ml 5% PBS-FBS 1400 rpm, 6 min  
Add ColdStain mix, incubate @4C for 30min  
Add 300-500 ul 1% RBC Lysis for 3 minutes  
Wash 2 ml 5% PBS-FBS 1400 rpm, 6min  
300 ul BD FixPerm, incubate @4C for 20min  
(vortex every 10 minutes)  
First PermWash: 1 ml PermWash 1500 rpm 6 min  
Second PermWash: 1 ml PermWash 1500 rpm 6 min  
Add Intracellular Stain, incubate @RT for 45min  
First PermWash: 2 ml PermWash 1500 rpm 6 min  
Resuspend in 70 ul 0.4% PFA-PBS  
Cap tubes, wrap rack in foil, store at 4C





## 3/7/2023

[illegible]



March 7th, 2023

Check for laser issues!


8:41: Template successfully set up!!

W15/V16  
V51/V16  
R6-R8

Fails QC earlier

start @ 8:47 pm

CD62L 1:07:50 min & 250k cells  
CD8 144k

Cord:   
shape

CD69 marks 150k events  
(this was ND050?)

//CCR4 stiller visible (not brightest)  
A lot of Vds ND050 !!  
decent cx on 3 stim

161 staining bit mean... adult.  
only slightly better cord

Δ intensity/shape to CD45RA cord/adult.  
... pure ATo.  
strong #NG

only 500 events for Vα24Ja18 (my suspicion why that first experiment failed so horribly)  
good thing remembered F10 CD3

Solid PD1 expression

strong but faint CD16 signal.

Faint Va7.2 signal

2 Zombie Cordobovs ... unstained signature? wrap it in.

Sc's wrapped @ 9:59 pm

1:10:00 minutes

2023\_ILT\_01

Samples @ 1001 pm

Info 18 0 00

1,200,000 cells Info 18 PMA 1,000,000  
1,300,000 cells Info 161 (solid monocytes)  
Info 180 PMA

Last sample  
10:38 pm

Revised 23%

Probably screen for dead-ish Inf  
@ Guava & select for 161  
Zombie?

minimal contains

Sc's 1 min

Samples 2 min

Info 18 0 00

|                      |        |         |
|----------------------|--------|---------|
| 3 million cells/ml   | 1 ml   | 2 fl    |
| 2 million cells/ml   | 666 fl | 1.33 fl |
| 1 million cells/ml   | 333 fl | 0.66 fl |
| 500k cells           | 167 fl | 0.33 fl |
| 2.5 million cells/ml | 833 fl | 1.66 fl |
| 1.5 million cells/ml | 500 fl | 1 fl    |

Samples - unstained @ 10:41 pm  
~400k cells (~330k)

ND006 - ctrl - SC - unstained - intracellular

ND006 - PMA - SC - unstained - intracellular

ND050 - ctrl - SC - unstained - intracellular

Adult unstained ctrl sc's no intracellular  
PMA

Shutdown @ 11:14 pm

# Unmixing Ab Experiments

★ NKG2D CD56 IFN $\gamma$  TNF $\alpha$  CD137 $\alpha$  ★



March 7th, 2023

2023\_ILT\_01

Initial Unmixing: unclear which cords used for what?

3 Dead → 018 } ctrl  
4 ~~Dead~~ Card ctrl → 101 }



March 7th, 2023

2023\_ILT\_01

