

November 25th, 2022

Intragen  
Nuclear  
FoxP3 Staining vs BD  
Fix Perm

Specimen	Status	Location	Conc	Date	Tasks	Volume	Ly	Ly+Mon	Total	2E+6	.3E+6
ND006	Adult		15E+6	03-04-21		2	5.35	6.71	10.7		
ND006	Adult		15E+6	03-04-21		2	6.13	7.59	12.26		
									23.0m		
						4	5.74	7.15		0.3484	
						5	4.60	5.72		0.434	65pl

Thaw spin @ 8:39 am  $\approx$  1300 rpm 10 min  
count stain @ 8:59 am  $\rightarrow$  9:06  
count done @ 9:20 am

Incubation @ 37°C start @ 9:52 am  $\rightarrow$  3:52 pm

Ab prep @ 2:24 pm

Done @ 3:19 pm

3:56 pm cells out.

4:29 aliquoted 2 in 1300/8 spin

15:46 pm  $\rightarrow$  cells into fridge  
not ideal, but lot sc's to deal with

Hi4a sc start 7:51 pm

(CCRB FoxP3) second up, had to take ~50pl  
from unstained, hindsight should've taken from  
V81 Fetc.

Also realized intracellular has Fetc FoxP3 so yeah.

Planning to spike in 0.5pl for the markers just in  
case they do matter...

Hot sc's @ 5:28 pm  $\rightarrow$  5:58 pm

Cold sc's @ 5:29 pm  $\rightarrow$  5:59 pm

Samples & LD sc @ 5:32 pm  $\rightarrow$  5:47 pm  
Spin @ 5:49 pm

2 sc's out @ 5:51 pm

53/54 Sp SB


59 00 63

Samples back @ 6:00

Ran at FBS

For stop, Ebtoscience mix @ 1pl ~ 500pl  
but instead of PMA  $\rightarrow$  substitute SEB (5pl/ml)  
for the activation condition.

Go for 6 hours




29 colors  
 $\times 2$   
58 sc's  
 $\times .3$   
17.4

23.0M  
17.4M  
5.60M

6 x 3 = 18  
52  
 $\times .3$   
15.6  $\rightarrow$  (16M)

2M SEB  
SCS



2M/500pl?

434pl + 66pl R10 = 500pl

2.5pl SEB + 1pl stop

2 tubes w/ CD107A samples

2 tubes activated SEB sc's + unstained

8 tubes ctrl tubes?

FBS: 68mls  
24mls  
132  
132  
132  
132

NucleoPerm wash 136mls

32  
 $\times 2$   
64  
136  
136  
136

2 samples

6mls FBS-PBS  
8mls Wash (4)  
1ml Fix nuclear

68  
 $\times 2$   
136  
136  
136

5:53 6:07  
6:16 pm sc's in FBS wash post RBC lysis

Samples into 37°C stain @ 6:20 pm → 6:50 pm

6:35 streptavidin sc's → 6:50 pm

sc's → stained ✓

lysed → ✓

(streptavidin ✓)

4°C @ 7:01 → 7:31

Fixed/Perm @ 7:19 → 7:20 → 7:30

6:37 Fix - 6:40 → 50 → 7:00 ✓

6:43 → 53 → 03 → 13 ✓

Perm Wash spin 1 ✓ 7:07 pm → 7:13 pm

2nd wash @ \_\_\_\_\_

nuc wash 2nd @ \_\_\_\_\_

7:31 pm

Rbc lysed @ 7:37

Just CO161 needs its final nuclear wash  
samples need streptavidin.

Samples streptavidin @ 7:47 pm → 8:02 pm

Sample streptavidin spin @ 8:04 pm

SC intracellular @ 8:18 pm <sup>+40</sup> → 8:58 pm

Two also mg 8:18

Sample Fixes @ 8:21 pm → 31 → 41 (FixPerm)

sc's 1st wash post stain @ 9:00 → 51 (Nuclear)  
last nuclear wash pre-intra for samples

Samples intra @ 9:09 → 9:49 pm

Final SC spin @ 9:15 pm

sc's away @ 9:34 pm

1st wash @ 9:50 pm ✓

2nd wash @ 9:58 pm ✓

Done @ 10:01



Unmixing ctrl																				
#	Filter	Single color (ul)	Ref ctrl type	name	Fluorochrome	Marker	Clone	Val lot #	During stim III	2	U/D RT 15 min	RT for 10 min	37°C for 30min	2	Surface 40C for 30min	2	4°C for 15 min	RG type then Nuclear/pe em	Inten Mix RT for 40min	2
1	UV2				BUV395	CD62L	(DREG-56)								1.2	2.4				
2	UV7				AF	AF-UV6														
3	UV9				BUV496	CD8	(RPA-18)								0.7	1.4			0.25	0.5
4	UV10				BUV563	CD69	(FN50)								0.5	1			0.5	1
5	UV11				BUV615	CCR4	(1G1)					2.5	5						0.5	1
6	U14				BUV661	Vd2	(B6)								0.7	1.4				
7	UV16				BUV737	CXCR3	(1C6/CXCR3)					2.5	5						0.5	1
8	V1				BUV805	CD4	(S3)								1.5	3			0.25	0.5
9	V3				BUV421	CD127	(A019P5)					1.5	3						0.5	1
10	V5				Pacific Blue	CD14	(M5E2)								2.0	4				
11	V7				Pacific Blue	CD19	(H1B19)								2.0	4				
12	V10				BV480	CD161-	(REA631)													
13	V11				BV570	Biotin														
14	V13				BV510	CD45RA	(H1100)													
15	V14				BV605	CXCR5	(U25204)					2.5	5						0.5	1
16	V15				BV650	CCR7	(G043H7)					1.5	3						0.5	1
17	B3				BV711	Va7.2	(3C10)								1.2	2.4			0.25	0.5
18	B4				BV750	IFN $\gamma$	(B27)												1.5	3
19	B6				BV786	CCR6	(11A9)					1.5	3							
20	B8				FOXP3/VD1	CD3	(SK7)													
21	B10				Spark Blue 550	PD-1	(EH12.2H7)								1.2	2.4			4	8
22	B13				PE	CD26	(M-A261)								1.2	2.4			0.25	0.5
23	B10				PE-Cy5	CD25	(M-A251)					1.2	2.4						0.25	0.5
24	B13				PerCP-Cy5.5	TNFA	(MAB11)								1.2	2.4			2.5	5
25	R1				PE-w0770	HLA-DR	(1243)													
26	R2				APC	CD39	(A1)													
27	R4				Alexa Fluor 647	IL2	(MO1-					1.2	2.4						1.2	2.4
28	R6				APC-R700	CD107a	(H4A3)													
29	R7				Zombie NIR	L/D														
30	R8				APC/Fire 750	CD27	(O323)					2	4							
31					APC/Fire 810	CD38	(HIT2)					1.6	3.2							
And UNMIXED CONTROLS !!!							Antibody Total	0.0			Antibody Total	18	36		14.9	29.8			14.7	29.4
							RT0 Media	20.5			Brilliant Stain	50	100		50	100			50	100
							Pipette draw volume /sample	19.5			Pipette draw volume /sample	64			60.9				60.7	

Notes: Need to resolve the what gets 37°C, what gets 4°C and what gets intracellular on Monday, can't be spiking two antibodies if avoidable for CD161.....

### Simplified Protocol

<Rest cells overnight in 6 well plate, 1E6/ml of R10>

**Collect, count, aliquot cells 2.0E+6 Cells/ 5ml polystyrene tube**

Cap and store incubator at 37°C for 2 two hours

**Add golgi block, and incubate for <6> hours**

Wash with 2 ml PBS, 1300rpm, 8min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6min

Add cold Surface mix, incubate @ 4C for 30min

Wash 2 ml 5% PBS-FBS 1500 rpm, 6min

1 ml Nuclear FixPerm, incubate @ 4C for 30 min

**First Perm Wash:** 2 ml NuclearWash 1500 rpm 6 min

**Second Perm Wash:** 2 ml NuclearWash 1500 rpm 6 min

**Add Intracellular Stain, incubate @ RT for 40min**

First Perm Wash: 2 ml NuclearWash 1500 rpm 6 min

Second Perm Wash: 2 ml NuclearWash 1500 rpm 6 min

Resuspend in 100  $\mu$ l 0.4% PFA-PBS

→ Not the LSR-II  
as reason ↓ 3pl  
for Harmon

Milky Biosesal  
PEAE31 or 14P-3610

or BD  
21  
Dceg5k

Testing on  
single cell  
data (Pik)

Supererogatory or not?



11/25/2022

[illegible]



November 25th, 2022

FoxP3 Staining

# Universal Intracellular Mix

eBioscience FoxP3/TF staining buffer set

00-5523

(Invitrogen/ThermoFisher)

FoxPerm Concentrate 4x 1:3 diluent

Buffer 10x w/ H<sub>2</sub>O

100  $\mu$ l Perm Buffer

x 2 ml wash

1 ml FoxP3 fixPerm 30 minutes

2 ml 1x  
x 2

(400-600 g 5 min)

2 ml 1x

Acquisition - Aurora on @ 9:58 am ✓

(moving my 2 Hastings files to backup disk / limited memory on Data :D)

Bailey is "using it for just an hour" sure, glances @ waste tank that I left full earlier that afternoon.

★ Running 50  $\mu$ l H<sub>2</sub>O @ 30 min mark, no cells in line.

60  
x 2 min  
120

QC @ 10:36 am (Lot 2004) Important as keeps gains set ~ spot/daily.

FixPerm start @ 10:41 am.

Unstaineds ← at whatever time frame they were at collecting 50  $\mu$ l / sec → ✓ for this round. (or 100,000 events /

CO<sub>2</sub> ↑ then used.

(CD69  
SEB)

0 ←  
more distinct doublets


11:00 am (stop to fix sampler ref names, Pump in correct collection amounts for samples that got extra cells)

(will need to just duplicate raw data to not need recreate again for FoxP3)

★ (don't think we have enough VdZ's) ★

... Did my SEB work? ... (did you vortex it?) ... IFN $\gamma$  no well ...  
if any very little?

From CD45 RA looks like we are having some major AF at least





November 25th, 2022

SEB that failed, or was it simultaneous addition of step?

FoxP3 Staining (did add step)

Va 7.2 similar, lost in the bell?

TNFA also looks  $\neq$ , what happened SEB? did I add these abs? I think so... culture condition same... and the marker times... slight...

or stop @ start?


If brighter signal is on a  $\Delta$  cell pop, does it matter if you focus exclusively on your own?

Lymphocytes are not no AF they have it, just swamped by monocyte levels. especially after culture? Y/N

Clear no cytokines or  $\Delta$  between tubes

~ some degradation ~

Particularly obvious single stain



healthy  
lost membrane integrity

1st CD62L bts below margin events, (nope) (some quality)

First set done @ 10:40  
12:19  
(1:40 min)


FoxP3 @ 12:25

Wow... vastly different. Almost at the threshold (250,000)

Lymphocyte average  $\sim$  1M multi

CD8 brighter w/ the FoxP3  $\sim$  higher AF background?

[set sds at 50  $\mu$ l volume]



but how does equivalent look compressed/spread over a log?

Not as easy as first glance

If you collect few events vs more events, effect on the final "color" signature sent to the matrix? [distinct lack of creativity gets on my nerves] ... I am side-whiplashing though...

What does a decrease in FSC actually reflect?

a decrease in internal cellular area from FoxP3? loss cytosol

[certain threshold % to be visualized on spectral line?]

They all picked up CD107a stain from culture period, so yeah, analogous to zombie NIR effectively



November 25th, 2022

FoxP3 Staining

If we do go nuclear stain route, we are going to need to change threshold, barely can gate on dead cells for its Zombie NIR SC.

(Swapping out waste before sample acquisition;

Sample acq @ 1:50 pm → Fixed 1.5 M FSC-A.

Manipulating FSC gain to have pop in same spot, but obscuring  $\Delta$  in properties


Do you want to spend rest of your life counting deer year after year - want to be up all night doing flow running cells? ... no.  
Leverage both to find useful things.

Done @ 2:09 → cleaning up

Unmixing 11/28/2022

FoxP3 → used so unstained for both act/unact unstained to start.

Last CD38 expression



No CCR3 for FixPerm

127 intrad

Matrix change of [0.05] as green cutoff.

Tag - on what?

48 62 6254.55  
-17 -22 2103  
31  
40 40.29

0.05 - 0.10 green //  
0.10 - 0.15 yellow //  
0.15 - 0.20 orange //  
0.20 - 0.30 red  
~~0.20 - 0.30~~  
~~0.30 - 0.35~~  
0.30 - 0.50 black





November 25th , 2022

FoxP3 Staining

Alexis meeting

Tuesday AM

Wednesday AM

Thursday AM

10:15 AM  
6:45  
17:00  
5pm  
Meeting 6pm

Tuesday  
5pm  
6pm

Wednesday  
5pm  
6pm

3-4pm  
Tuesday 10pm  
Tuesday 10pm

- Kirschen meeting
- Mancini/Andrew
- Alexis Gonzales