

August 23, 2021

HEU: Innate-like T cells
Mixing Tetramers & Biotin

Specimen	Status	Location	Conc	Date	Tasks	Volume	Ly	Ly+M	Total	1.E+6	2E+6
ND006						3	3.78E4			265pl	
ND006						3	3.93E4			254pl	
CQ1070 #3611						3	5.3E4			189pl	

[6 tubes]

Tetramer Test	Ex vivo Both	Surface μL	x 6.5
V450 BV421	CD4	1.8	11.7
V525 BV510	Aqua L/D	1:500	-
V670 BV650			
B530 Alexa488 (FITC)	hCD1d Tetramer		
B710 PerCPeF710	CD45RA	2	13
Y590 PE	CD161	0.5	3.3
Y615 PE Dazzle			
Y780 PE Vio770	CD3	0.5	3.3
R670 APC	MR1 Tetramer		
R780 APCFire750	CD8	1.5	9.75
20.5 $\mu\text{L}/\text{rxn}$	PBS	11.2	72.8

(17.5 draw)

← 2pl
← leaving these at master mix

← 5pl

2 draws each Fite hCD1d/solo
12 each APC MR1/CD45RA/CD8
x 5
60pl
6pl tet into 54 (1:10)

22M
Pib cont Pib's biot cont cont
Aqua 1 6 - 6 2 4 = 19
Horizon 3 - 6 - - - = 9
x 200 = 15,200
9pl in 9mls
30.4pl in 15.2mls
32
x 3
166

LOR messed up CD3 FSC gate (CD misread)

Tetramer Test	Ex vivo Both	Surface μL	x
V450 BV421	CD4	1.8	
V525 BV510	Aqua L/D	1:500	
V670 BV650			
B530 Alexa488 (FITC)	hCD1d Tetramer		
B710 PerCPeF710	CD45RA	2	
Y590 PE	CD161	0.5	
Y615 PE Dazzle			
Y780 PE Vio770	CD3	0.5	
R670 APC	MR1 Tetramer		
R780 APCFire750	CD8	1.5	
20.5 $\mu\text{L}/\text{rxn}$	PBS		

Tetramer Test	Ex vivo Both	Surface μL	x
V450 BV421	CD4	1.8	
V525 BV510	Aqua L/D	1:500	
V670 BV650			
B530 Alexa488 (FITC)	hCD1d Tetramer		
B710 PerCPeF710	CD45RA	2	
Y590 PE	CD161	0.5	
Y615 PE Dazzle			
Y780 PE Vio770	CD3	0.5	
R670 APC	MR1 Tetramer		
R780 APCFire750	CD8	1.5	
20.5 $\mu\text{L}/\text{rxn}$	PBS		

Done @ 4:50pm

August 23, 2021

HEU: Innate-like T cells
Mixing Tetramers & Biotin NKT P2 + CCR7
w/ NKG2D ~ biotin non biotin

6 tubes?}

master mix

Vα Test A	Ex vivo NKT-p2	Surface μL x 1	Intra μL x 4.5
BV421	Va24	2	
BV510	CD3	1.5	
BV650	NKG2D	2 (1.5)	
	Streptavidin		
Alexa 488 (FITC)	dg9	-	3 13.5
PerCPeF710	CD56	1	
PE	PD1	1.5	
PE Dazzle			
PE Vio770	NKG2A	0.6	
APC	CCR6	2.5	
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$; PBS		9.4	18.5 83.3

Vα Test B	Ex vivo NKT-p2	Surface μL x 1	Intra μL x
BV421	Va24	2	
BV510	CD3	1.5	
BV650	NKG2D	2 (1.5)	
	Streptavidin		
Alexa 488 (FITC)	dg9	-	3
PerCPeF710	CD56	1	
PE	PD1	1.5	
PE Dazzle	CCR7	1.5	
PE Vio770	NKG2A	0.6	
APC	CCR6	2.5	
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$; PBS		7.9	

* I added right?
18.5
x 4.5
83.3
74.0
83.25

Tetramer A	Ex vivo NKT-p2	Surface μL x 1	Intra μL x
BV421	Va24	2	
BV510	CD3	1.5	
BV650	NKG2D? fixed	1.0 \downarrow	
Alexa 488 (FITC)	dg9	-	3
PerCPeF710	CD56	1	
PE	PD1	1.5	
PE Dazzle	CCR7	1.5	
PE Vio770	NKG2A	0.6	
APC	hCD1d PBS-57/solo		
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$; PBS		14.9	

Tetramer B	Ex vivo NKT-p2	Surface μL x 1	Intra μL x
BV421	Va24	1.5	
BV510	CD3	1.5	
BV650	NKG2D? fixed	1.0 \downarrow	
Alexa 488 (FITC)	dg9	-	3
PerCPeF710	CD56	1	
PE	PD1	1.5	
PE Dazzle	CCR7	1.5	
PE Vio770	NKG2A	0.6	
APC	hCD1d PBS-57/solo		
APC Fire750	Horizon L/D	1:1000	
20.5 $\mu\text{L}/\text{rxn}$; PBS		13.4	

RBC lysis

Sheet 2 aliquots

Can't remember if I grabbed biotin NKG2D or CD26... I think safe? well find out...

* Healthy Vα staining vs CCR6 in this panel.

* FSC off because you fixed primed these, switched 59000 lot dead cells in solo A tet

Everything is messy as hell.

August , 2021

HEU: Innate-like T cells Mixing Tetramers & Biotin


- 1) Tetramer 40 min RT, wash, Abs 20 min 4°C wash, Streptavidin 15 minutes 4°C, Wash
- 2) Tetramer + Abs 20 min 4°C, Strep 15 min [known horror]
- 3) CD26 + Strep, Wash, Tetramer 40 min RT, wash, Abs 20 min 4°C wash

Tetramer Test B1	Ex vivo APC	Surface μL x 11	
V450 BV421	CCR6	20	27.5
V525 BV510	Aqua L/D	1:50	
V670 BV650	CD26 - Strep	0	
B530 FITC	CD4	1.8	19.8
B710 PerCPeF710	CD45RA	2	22
Y590 PE	CD161	0.5	5.5
Y615 PE Dazzle	CCR7	1.5	16.5
Y780 PE Vio770	CD3	0.5	5.5
R670 APC	Tetramer (+/-)		
R780 APCFire750	CD8	1.5	16.5
20.5 $\mu\text{L}/\text{rxn}$	PBS	5.2	57.2

could see if ~~that~~ original way maybe working too.

→ 22 (5.5 pB3)

yes was done
leaving this at max time w/ APC-stain



May not be as distinguishable since MAFs are robust.

(5) 15.5

~~15.5 draw~~

(15.5 draw)

(15.5 draw)

2 ml washes in between Strep containing steps

$$\begin{array}{r} 300 \\ \times 8 \\ \hline 3000 \end{array}$$

P200s out intra @ 14pm

Tailing step starts out @ same time

P2 Va VB into intra @ 3:48 → 4:22 pm out → Done

P2 ~~for~~ AB into first spin

Same wash @ 1400/6

Remaining P2 in @ 4:06 → 4:46 done

August 23, 2021

HEU: Innate-like T cells
Mixing Tetramers & Biotin

APC MALT [20 min 4°C] vs [40 min RT + 20 min 4°C] {4 tubes}

Tetramer Test B1	Ex vivo APC	Surface μ L x	
V450 BV421	CCR6	2.5	
V525 BV510	Aqua L/D	1:500	
V670 BV650			
B530 FITC	CD4	1.8	
B710 PerCPeF710	CD45RA	2	
Y590 PE	CD161	0.5	
Y615 PE Dazzle	CCR7	1.5	
Y780 PE Vio770	CD3	0.5	
R670 APC	Tetramer		
R780 APCFire750	CD8	1.5	
20.5 μ L/rxn;	PBS		

[postpaul]

* APC OP-RV / GFP in cord

* APC hCD1d / solo in cord

- baseline characterization. (abs ✓ just more tetramer)

order more CD3 BULG
* CCR6 BV421 = APC caps
Very similar !!

LSR 08/24/21 → Check to setup a run Dazzle comp

NOOK MALTs
Centrals ~ 1% CD3 CCR6 + CD45RA - CCR7 - CD8
NK1T staining ← CD161, real, wrong marker? TBD...

OP-RV CQ vly OP-RV ^{smear} on CD3 - pop ... my findings → BV650? nope on APC
after collection graph revealed.

none on hCD1d, visible population.... mistaken tetramers?
just compensation is shifted. Few CD4s.

2.21	1/a	2
50	1/D	Aqua
10	-	
10	CO4	1.8
10	CO5DA	2
10	CO161	0.5
10	CCR7	1.5
10	CD3	0.5
10	CCR6	2.5
10	CO8	1.5
10	PBS	8.2

Single Va 24

Va 7.2
mix controls

~~add in CO2...~~

after surface
strain ...
created
decision